

COMPLEX III SUBUNIT 7 RABBIT MAB

货号: N263055

产品全名: Complex III Subunit 7 兔单克隆抗体

基因符号 QPC; QCR7; QP-C; UQBC; UQBP; UQPC; UQCR6; MC3DN3

UNIPROT ID: P14927

背景: This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This component is involved in redox-linked proton pumping. Caution Was originally thought to be the ubiquinone-binding protein (QP-C).

抗原: Recombinant protein of human UQCRB

经过测试的应用: WB, IHC-P, IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-4F5

分子量: Calculated MW: 14 kDa; Observed MW: 14 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using UQCRB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of UQCRB in rat Brain lysates using Complex III Subunit 7 antibody.