

CRABP2 RABBIT MAB

货号: N261588

产品全名: CRABP2 兔单克隆抗体

基因符号 Cellular retinoic acid-binding protein II; CRABP-II

UNIPROT ID: P29373

背景: Transports retinoic acid to the nucleus.

抗原: Recombinant protein of human CRABP2

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-4C4

分子量: Calculated MW: 16 kDa; Observed MW: 16 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of CRABP2 (green) in MCF-7 using CRABP2 antibody, and DAPI (blue)



Western blot analysis of CRABP2 in HeLa lysates using CRABP2 antibody.