

CRIPTO 3 RABBIT PAB

货号: N225205

产品全名: Cripto 3 兔多抗

基因符号 TDGFIP3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3; Cripto-3 growth factor; Epidermal growth factor-like cripto protein CR3; Teratocarcinoma-derived growth factor 1 pseudogene 3

UNIPROT ID: P51864

背景: Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human TDGFIP3. AA range:1-50

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 21 kDa; Observed MW: 21 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cripto 3 in HepG2, 293T, AD293, HeLa lysates using Cripto 3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
