

CTGF RABBIT MAB

货号: N262074

产品全名: CTGF 兔单克隆抗体

基因符号 CCN2; NOV2; HCS24; IGFBP8

UNIPROT ID: P29279

背景: Major connective tissue mitogen secreted by vascular endothelial cells. Promotes proliferation and differentiation of chondrocytes. Mediates heparin- and divalent cation-dependent cell adhesion in many cell types including fibroblasts, myofibroblasts, endothelial and epithelial cells. Enhances fibroblast growth factor-induced DNA synthesis.

抗原: Recombinant protein of human CTGF

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-8C2

分子量: Calculated MW: 38 kDa; Observed MW: 35 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CTGF in C6, HeLa lysates using CTGF antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
