

CULLIN 4A RABBIT PAB

货号: N225946

产品全名: Cullin 4A 兔多抗

基因符号 Cul4a; Cullin-4A

UNIPROT ID: Q13619

背景: Core component of multiple cullin-RING-based E3 ubiquitin-protein ligase complexes which mediate the ubiquitination and subsequent proteasomal degradation of target proteins. As a scaffold protein may contribute to catalysis through positioning of the substrate and the ubiquitin-conjugating enzyme.

抗原: A synthesized peptide derived from human Cullin4A

经过测试的应用: WB,IHC-P,ICC/IF,IP,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 88 kDa; Observed MW: 88 kDa

亚型: IgG

纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cullin 4a in HeLa lysates using Cullin 4A antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
