

CYTOCHROME P450 REDUCTASE RABBIT PAB

货号: N225965

产品全名: Cytochrome P450 Reductase 兔多抗

基因符号 CPR; CYPOR; P450R; por

UNIPROT ID: P16435

背景: This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5.

抗原: A synthesized peptide derived from human Cytochrome P450 Reductase

经过测试的应用: WB,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 77 kDa; Observed MW: 77 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cytochrome P450 Reductase in A431 lysates using Cytochrome P450 Reductase antibody.