

## CYTOKERATIN 9 RABBIT PAB

货号: N225713

产品全名: Cytokeratin 9 兔多抗

基因符号 K9; CK-9; EPPK

**UNIPROT ID:** P35527

**背景:** May serve an important special function either in the mature palmar and plantar skin tissue or in the morphogenetic program of the formation of these tissues. Plays a role in keratin filament assembly.

**抗原:** A synthesized peptide derived from human Cytokeratin 9

**经过测试的应用:** WB,IHC-P,IP

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 62 kDa; Observed MW: 62 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human, Mouse and Rat

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cell Biology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cytokeratin 9 in HeLa lysates using Cytokeratin 9 antibody.