

## DAPI1 RABBIT MAB

货号: N262112

产品全名: DAPI1 兔单克隆抗体

基因符号 Death-associated protein 1; DAP-1

**UNIPROT ID:** P51397

背景: Negative regulator of autophagy. Involved in mediating interferon-gamma-induced cell death.

抗原: A synthetic peptide of human DAPI1

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-2J5

分子量: Calculated MW: 11 kDa; Observed MW: 11 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DAPI1 in MCF-7, C6, 3T3 lysates using DAPI1 antibody.