

DDX3 MOUSE MAB

货号: N261537

产品全名: DDX3 小鼠单抗

基因符号 DBX; DDX3; HLP2; DDX14.

UNIPROT ID: O00571

背景: Multifunctional ATP-dependent RNA helicase.

抗原: Purified recombinant human DDX3 protein fragments expressed in E.coli.

经过测试的应用: WB,IP,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 6G8-F4-E3

分子量: Calculated MW: 73 kDa; Observed MW: 73 kDa

亚型: IgG2a


纯化: Affinity Purified

种属反应性: Human,Rat,Mouse,Monkey

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DDX3 in COS7, Jurkat, C6, 3T3 and HeLa lysates using DDX3 antibody.