

DEK RABBIT MAB

货号: N262125

产品全名: DEK 兔单克隆抗体

基因符号 Dek; DEK oncogene; Protein DEK

UNIPROT ID: P35659

背景: Involved in chromatin organization.

抗原: Recombinant protein of human DEK

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-5D2

分子量: Calculated MW: 43 kDa; Observed MW: 50 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using DEK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of DEK in K562, MCF-7, HeLa lysates using DEK antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
