

DESMOCOLLIN 3 RABBIT MAB

货号: N261610

产品全名: Desmocollin 3 兔单克隆抗体

基因符号 Cadherin family member 3; Desmocollin-4; HT-CP

UNIPROT ID: Q14574

背景: Component of intercellular desmosome junctions.

抗原: A synthetic peptide of human Desmocollin 3

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-4G2

分子量: Calculated MW: 100 kDa; Observed MW: 100 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Desmocollin 3 in K562, HeLa lysates using Desmocollin 3 antibody.



Immunocytochemistry analysis of Desmocollin 3 (green) in HeLa using Desmocollin 3 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
