

## DLX4 RABBIT MAB

货号: N262139

产品全名: DLX4 兔单克隆抗体

基因符号 DLX4; BPI; DLX7; DLX8; DLX9; Homeobox protein DLX-4; Beta protein 1; Homeobox protein DLX-7; Homeobox protein DLX-8

**UNIPROT ID:** Q92988

背景: May play a role in determining the production of hemoglobin S. May act as a repressor. During embryonic development, plays a role in palatogenesis.

抗原: A synthetic peptide of human Dlx4

经过测试的应用: WB,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-1C5

分子量: Calculated MW: 26 kDa; Observed MW: 35 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Dlx4 in K562, rat Brain, HeLa lysates using Dlx4 antibody.



Immunocytochemistry analysis of DLX4 (green) in K562 using DLX4 antibody, and DAPI (blue).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---