

DNAJA4 RABBIT MAB

货号: N262141

产品全名: DNAJA4 兔单克隆抗体

基因符号 MST104; MSTP104; PRO1472

UNIPROT ID: Q8WW22

背景: Enables chaperone binding activity and unfolded protein binding activity. Involved in several processes, including negative regulation of endothelial cell migration; negative regulation of inclusion body assembly; and protein refolding. Located in cytosol and membrane.

抗原: A synthetic peptide of human DNAJA4

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-8G5

分子量: Calculated MW: 45 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DNAJA4 in C6, rat Brain lysates using DNAJA4 antibody.



HeLa lysates using DNAJA4 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
