

## DNAJC8 RABBIT MAB

货号: N262136

产品全名: DNAJC8 兔单克隆抗体

基因符号 SPF31; HSPC331

**UNIPROT ID:** O75937

背景: Suppresses polyglutamine (polyQ) aggregation of ATXN3 in neuronal cells (PubMed:27133716).

抗原: A synthetic peptide of human DjC8

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-311

分子量: Calculated MW: 30 kDa; Observed MW: 30 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat,Monkey


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DjC8 in Jurkat, Hela, COS7, Raji lysates using DjC8 antibody.



Western blot analysis of DjC8 in mouse brain, rat lung lysates using DNAJC8 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---