

EIF4A1 RABBIT MAB

货号: N262167

产品全名: eIF4A1 兔单克隆抗体

基因符号 Eukaryotic initiation factor 4A-I; eIF-4A-I; eIF4A-I; ATP-dependent RNA helicase eIF4A-1

UNIPROT ID: P60842

背景: ATP-dependent RNA helicase which is a subunit of the eIF4F complex involved in cap recognition and is required for mRNA binding to ribosome.

抗原: A synthetic peptide of human eIF4A1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-2K6

分子量: Calculated MW: 46 kDa; Observed MW: 46 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of eIF4A1 in K562, C6, 3T3, HeLa lysates using eIF4A1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
