

EMERIN RABBIT PAB

货号: N225870

产品全名: Emerin 兔多抗

基因符号 STA; EDMD; LEMD5

UNIPROT ID: P50402

背景: Emerin may regulate gene expression through binding to other transcriptional regulators. Emerin binds to β -catenin and inhibits its nuclear accumulation. Recent studies demonstrate that Emerin is required for HIV-1 infectivity.

抗原: A synthesized peptide derived from human Emerin

经过测试的应用: WB,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 29 kDa; Observed MW: 29 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing



Western blot analysis of Emerin in 293T lysates using Emerin antibody.