

EYA2 RABBIT MAB

货号: N263323

产品全名: EYA2 兔单克隆抗体

基因符号 EAB1

UNIPROT ID: O00167

背景: Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.

抗原: A synthetic peptide of human EYA2

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-6A1

分子量: Calculated MW: 59 kDa; Observed MW: 59 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.