

F4/80 RAT MAB

货号: N263636

产品全名: F4/80 Rat 单抗

基因符号 Emrl; Ly71; Gpf480; TM7LN3; DD7A5-7; EGF-TM7

UNIPROT ID: Q6I549

背景: Orphan receptor involved in cell adhesion and probably in cell-cell interactions specifically involving cells of the immune system. May play a role in regulatory T-cells (Treg) development.

抗原: A synthesized peptide derived from Mouse F4/80.

经过测试的应用: IHC-P, ICC/IF, FC

推荐稀释比: IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rat

克隆性: Rat Monoclonal

克隆编号: 31B1

分子量: -

亚型: IgG2a

纯化: Affinity Chromatography

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Microglia

储存和运输: Store at -20°C. Avoid repeated freezing and thawing