

FDFT1 RABBIT MAB

货号: N261639

产品全名: FDFT1 兔单克隆抗体

基因符号 DGPT; ERG9; FDFT1; SQS; Squalene synthase; SS

UNIPROT ID: P37268

背景: Critical branch point enzyme of isoprenoid biosynthesis that is thought to regulate the flux of isoprene intermediates through the sterol pathway.

抗原: Recombinant protein of human FDFT1

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-6B1

分子量: Calculated MW: 48 kDa; Observed MW: 48 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of FDFT1 (green) in HeLa using FDFT1 antibody, and DAPI (blue)



Western blot analysis of FDFT1 in K562, C6, HeLa lysates using FDFT1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
