

## FETUIN A RABBIT MAB

货号: N262216

产品全名: Fetuin A 兔单克隆抗体

基因符号 16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743

**UNIPROT ID:** P02765

背景: Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

抗原: A synthetic peptide of human Fetuin A

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-5A4

分子量: Calculated MW: 39 kDa; Observed MW: 45-50 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---