

FGF11 RABBIT PAB

货号: N225191

产品全名: FGF11 兔多抗

基因符号 FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3

UNIPROT ID: Q92914

背景: Probably involved in nervous system development and function.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human FGF11. AA range:21-70

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 25 kDa; Observed MW: 25 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.