

FGFR3 RABBIT MAB

货号: N263132

产品全名: FGFR3 兔单克隆抗体

基因符号 FGFR3; JTK4; Fibroblast growth factor receptor 3; FGFR-3; CD antigen CD333

UNIPROT ID: P22607

背景: Receptor for acidic and basic fibroblast growth factors. Preferentially binds FGF1.

抗原: A synthetic peptide of human FGFR3

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-6F4

分子量: Calculated MW: 88 kDa; Observed MW: 125 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of FGFR3 (green) in HeLa using FGFR3 antibody, and DAPI (blue)



Western blot analysis of FGFR3 in HeLa lysates using FGFR3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using FGFR3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.