

FIG4 RABBIT MAB

货号: N263127

产品全名: FIG4 兔单克隆抗体

基因符号 YVS; BTOP; SAC3; ALS11; CMT4J; KIAA0274; dJ24914.1

UNIPROT ID: Q92562

背景: Dual specificity phosphatase component of the PI(3,5)P2 regulatory complex which regulates both the synthesis and turnover of phosphatidylinositol 3,5-bisphosphate (PtdIns(3,5)P2).

抗原: A synthetic peptide of human FIG4

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-4C6

分子量: Calculated MW: 104 kDa; Observed MW: 104 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FIG4 in HL-60, U2OS, U87-MG lysates using FIG4 antibody.