

FKBP51 RABBIT MAB

货号: N262223

产品全名: FKBP51 兔单克隆抗体

基因符号 P54; AIG6; FKBP51; FKBP54; PPIase; Ptg-10

UNIPROT ID: Q13451

背景: Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

抗原: A synthetic peptide of human FKBP51

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-8E5

分子量: Calculated MW: 51 kDa; Observed MW: 51 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Hamster

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.