

## FKBP52 RABBIT MAB

货号: N262224

产品全名: FKBP52 兔单克隆抗体

基因符号 HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PPIase

**UNIPROT ID:** Q02790

背景: Immunophilin protein with PPIase and co-chaperone activities (By similarity). Component of unliganded steroid receptors heterocomplexes through interaction with heat-shock protein 90 (HSP90).

抗原: A synthetic peptide of human FKBP52

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-4A3

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat, Hamster


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of FKBP52 (green) in HeLa using FKBP52 antibody, and DAPI (blue)



Western blot analysis of FKBP52 in Raw264.7 lysates using FKBP52 antibody.



Western blot analysis of FKBP52 in Jurkat, C6, CHO-K1, HeLa lysates using FKBP52 antibody.