

FLOTILLIN 1 RABBIT MAB

货号: N262226

产品全名: Flotillin 1 兔单克隆抗体

基因符号 Flotillin-1

UNIPROT ID: O75955

背景: May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

抗原: A synthetic peptide of human Flotillin 1

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-5E9

分子量: Calculated MW: 47 kDa; Observed MW: 47 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Flotillin 1 (green) in HeLa using Flotillin 1 antibody, and DAPI (blue)



Western blot analysis of Flotillin 1 in K562, HeLa lysates using Flotillin 1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
