

## FLT3 LIGAND RABBIT PAB

货号: N225190

产品全名: Flt3 Ligand 兔多抗

基因符号 FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine

**UNIPROT ID:** P49771

背景: Stimulates the proliferation of early hematopoietic cells by activating FLT3.

抗原: The antiserum was produced against synthesized peptide derived from the C-terminal region of human FLT3LG. AA range:171-220

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 26 kDa; Observed MW: 26 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Flt3 Ligand in HeLa lysates using Flt3 Ligand antibody.