

FOXA1 RABBIT MAB

货号: N262232

产品全名: FOXA1 兔单克隆抗体

基因符号 FOXA1; Forkhead box A1; Forkhead box protein A1; HNF-3-alpha; TCF3A; HNF3A; Transcription factor 3A

UNIPROT ID: P55317

背景: Forkhead box protein A1 is a transcription factor required for the development of endoderm-derived organs, such as liver, lung, and prostate. FoxA1 functions as a pioneer factor that is recruited primarily to the distant enhancers to change chromatin structure for transcription in a cell type-specific manner. Involved in regulation of apoptosis by inhibiting the expression of BCL2. Involved in cell cycle regulation by activating expression of CDKN1B, alone or in conjunction with BRCA1.

抗原: Recombinant protein of human FOXA1

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-5A7

分子量: Calculated MW: 49 kDa; Observed MW: 49 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of FOXA1 (green) in MCF-7 using FOXA1 antibody, and DAPI (blue)



Western blot analysis of FOXA1 in T47D lysates using FOXA1 antibody.



Western blot analysis of FOXA1 in C6 lysates using FOXA1 antibody.