

FOXP1 RABBIT MAB

货号: N262236

产品全名: FOXP1 兔单克隆抗体

基因符号 FOXP1; HSPC215; Forkhead box protein P1

UNIPROT ID: Q9H334

背景: Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTBP1 is not essential. Essential transcriptional regulator of B cell development.

抗原: A synthetic peptide of human FOXP1

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-613

分子量: Calculated MW: 75 kDa; Observed MW: 82-90 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FOXP1 in Raw264.7, rat lung lysates using FOXP1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
