

FRAGILIS RABBIT PAB

货号: N225888

产品全名: Fragilis 兔多抗

基因符号 IFITM3; Interferon-induced transmembrane protein 3;
Dispanin subfamily A member 2b; DSPA2b; Interferon-inducible protein 1-8U

UNIPROT ID: Q01628

背景: IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens, namely influenza A H1N1 virus, West Nile virus (WNV), and dengue virus (WNV), by inhibiting the early step(s) of replication.

抗原: A synthesized peptide derived from human Fragilis

经过测试的应用: WB,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 15 kDa; Observed MW: 15 kDa

亚型: IgG

纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Stem Cells

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Fragilis in HeLa lysates using Fragilis antibody.