

FRS2 RABBIT PAB

货号: N225844

产品全名: FRS2 兔多抗

基因符号 Fibroblast growth factor receptor substrate 2; FGFR substrate 2; FGFR-signaling adaptor SNT; Sucl-associated neurotrophic factor target 1; SNT-1

UNIPROT ID: Q8WU20

背景: Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.

抗原: A synthesized peptide derived from human FRS2

经过测试的应用: WB,IP,FC

推荐稀释比: WB: 1/500-1/1000 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 57 kDa; Observed MW: 57 kDa

亚型: IgG

纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FRS2 in HeLa lysates using FRS2 antibody.