

FSP1 RABBIT PAB

货号: N225430

产品全名: FSP1 兔多抗

基因符号 AIFM2; AMID; PRG3; Apoptosis-inducing factor 2; Apoptosis-inducing factor homologous mitochondrion-associated inducer of death; Apoptosis-inducing factor-like mitochondrion-associated inducer of death; p53-responsive gene 3 蛋白

UNIPROT ID: Q9BRQ8

背景: A NAD(P)H-dependent oxidoreductase involved in cellular oxidative stress response.

抗原: Synthesized peptide derived from the Internal region of human AMID. at AA range: 110-190

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 41 kDa; Observed MW: 38 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human,Mouse,Monkey

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FSP1 in COS7 lysates using FSP1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
