

FUNDC1 RABBIT PAB

货号: N225130

产品全名: FUNDC1 兔多抗

基因符号 FUNDC1; FUN14 domain-containing protein 1

UNIPROT ID: Q8IVP5

背景: Acts as an activator of hypoxia-induced mitophagy, an important mechanism for mitochondrial quality control.

抗原: The antiserum was produced against synthesized peptide derived from the N-terminal region of human FUNDC1. AA range:1-50

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 17 kDa; Observed MW: 17 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FUNDC1 in AD293 lysates using FUNDC1 antibody.



Western blot analysis of FUNDC1 in AD293 lysates using FUNDC1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
