

GBP1 RABBIT MAB

货号: N263342

产品全名: GBP1 兔单克隆抗体

基因符号 GTP-binding protein 1

UNIPROT ID: P32455

背景: Hydrolyzes GTP to GMP in 2 consecutive cleavage reactions (By similarity).

抗原: A synthetic peptide of human GBP1

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-7B6

分子量: Calculated MW: 68 kDa; Observed MW: 68 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using GBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of GBP1 (green) in Hela using GBP1 antibody, and DAPI (blue).



Western blot analysis of GBP1 in HepG2, U251, HL-60, Lncap lysates using GBP1 antibody.