

GIV RABBIT PAB

货号: N226033

产品全名: GIV 兔多抗

基因符号 Akt phosphorylation enhancer; APE; Coiled-coil domain-containing protein 88A; G alpha-interacting vesicle-associated protein; GIV; HkRP1

UNIPROT ID: Q3V6T2

背景: Bifunctional modulator of guanine nucleotide-binding proteins (G proteins).

抗原: Recombinant protein of human Girdin

经过测试的应用: WB, ICC/IF, FC, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 216 kDa; Observed MW: 216 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Girdin in HeLa, A549, HL-60, U2OS, C6 lysates using Girdin antibody.



Western blot analysis of GIV in mouse brain, Raw264.7 lysates using GIV antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
