

GLI3 RABBIT MAB

货号: N262277

产品全名: GLI3 兔单克隆抗体

基因符号 GLI3-190; GLI3FL; GLI3-83; GLI3 C-terminally truncated form

UNIPROT ID: Q5IS56

背景: Has a dual function as a transcriptional activator and a repressor of the sonic hedgehog (Shh) pathway, and plays a role in limb development.

抗原: Recombinant protein of Gli3

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2B9

分子量: Calculated MW: 170 kDa; Observed MW: 200 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Gli3 in 293 lysates using GLI3 antibody.