

## GLUTATHIONE REDUCTASE RABBIT MAB

货号: N262284

产品全名: Glutathione Reductase 兔单克隆抗体

基因符号 GR; HEL-75; HEL-S-122m

**UNIPROT ID:** P00390

背景: Maintains high levels of reduced glutathione in the cytosol.

抗原: A synthetic peptide of human Glutathione Reductase

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-2E4

分子量: Calculated MW: 56 kDa; Observed MW: 51 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Glutathione Reductase in HeLa, Jurkat lysates using Glutathione Reductase antibody.