

GLUTATHIONE SYNTHETASE RABBIT MAB

货号: N262287

产品全名: Glutathione Synthetase 兔单克隆抗体

基因符号 Glutathione synthase; Glutathione synthetase; GSH S; GSH synthetase; GSH-S; GSHS; GSS; HEL-S-64p; HEL-S-88n

UNIPROT ID: P48637

背景: Sulfur metabolism; glutathione biosynthesis; glutathione from L-cysteine and L-glutamate: step 2/2.

抗原: A synthetic peptide of human Glutathione Synthetase

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-5G5

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Glutathione Synthetase in HeLa, mouse brain, 293T, Jurkat lysates using Glutathione Synthetase antibody.



Western blot analysis of Glutathione Synthetase in rat kidney, rat liver lysates using Glutathione Synthetase antibody.