

## GNE RABBIT MAB

货号: N262276

产品全名: GNE 兔单克隆抗体

基因符号 NM; DMRV; IBM2; Uae1; GLCNE

**UNIPROT ID:** Q9Y223

**背景:** Regulates and initiates biosynthesis of N-acetylneuraminic acid (NeuAc), a precursor of sialic acids. Plays an essential role in early development. Required for normal sialylation in hematopoietic cells. Sialylation is implicated in cell adhesion, signal transduction, tumorigenicity and metastatic behavior of malignant cells.

**抗原:** A synthetic peptide of human GLCNE

**经过测试的应用:** WB,IP

**推荐稀释比:** WB: 1/500-1/1000 IP: 1/20

**种属反应性:** Rabbit

**克隆性:** Rabbit Monoclonal

**克隆编号:** R05-7A8

**分子量:** Calculated MW: 79 kDa; Observed MW: 75 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Human

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Cardiovascular

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GNE in K562, HeLa lysates using GNE antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---