

GRK3 RABBIT PAB

货号: N225956

产品全名: GRK3 兔多抗

基因符号 ADRBK2; BARK2; GRK3

UNIPROT ID: P35626

背景: Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

抗原: A synthesized peptide derived from human GRK3

经过测试的应用: WB,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 80 kDa; Observed MW: 80 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GRK3 in Ramos lysates using GRK3 antibody.