

## GRK3 RABBIT PAB

货号: N226028

产品全名: GRK3 兔多抗

基因符号 ADRBK2; BARK2; GRK3

**UNIPROT ID:** P35626

背景: Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

抗原: A synthetic peptide of human GRK3

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 80 kDa; Observed MW: 80 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using GRK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GRK3 in Jurkat lysates using GRK3 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---