

## GROWTH HORMONE RECEPTOR RABBIT MAB

货号: N262303

产品全名: Growth Hormone Receptor 兔单克隆抗体

基因符号 GHR; Growth hormone receptor; GH receptor; Somatotropin receptor

**UNIPROT ID:** P10912

背景: Receptor for pituitary gland growth hormone involved in regulating postnatal body growth. On ligand binding, couples to the JAK2/STAT5 pathway. The soluble form (GHBP) acts as a reservoir of growth hormone in plasma and may be a modulator/inhibitor of GH signaling. Isoform 2 up-regulates the production of GHBP and acts as a negative inhibitor of GH signaling.

抗原: A synthetic peptide of human Growth hormone receptor

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-6E6

分子量: Calculated MW: 72 kDa; Observed MW: 120 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Growth hormone receptor in rat Brain lysates using Growth Hormone Receptor antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---