

GSDME RABBIT MAB

货号: N262128

产品全名: GSDME 兔单克隆抗体

基因符号 DFNA5; ICERE-1; GSDME

UNIPROT ID: O60443

背景: Plays a role in the TP53-regulated cellular response to DNA damage probably by cooperating with TP53 (PubMed:16897187, PubMed:18223688).

抗原: Recombinant protein of human DFNA5/GSDME

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-5II

分子量: Calculated MW: 55 kDa; Observed MW: 55 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DFNA5/GSDME in U251MG, U2OS, HepG2, HeLa lysates using GSDME antibody.