

## GST3 RABBIT MAB

货号: N262309

产品全名: GST3 兔单克隆抗体

基因符号 GSTP1; FAEES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1

**UNIPROT ID:** P09211

背景: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

抗原: A synthetic peptide of human GST3/GST pi

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-2F4

分子量: Calculated MW: 23 kDa; Observed MW: 23 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using GST3/GST pi antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GST3/GST pi in HepG2, Hela lysates using GST3 antibody.