

GST3 RABBIT PAB

货号: N225988

产品全名: GST3 兔多抗

基因符号 GSTP1; FAES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1

UNIPROT ID: P09211

背景: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

抗原: A synthesized peptide derived from human GST3 / GST pi

经过测试的应用: WB,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 23 kDa; Observed MW: 23 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GST3 / GST pi in PC-3 lysates using GST3 antibody.