

## **GTF2H1 (3E8) MOUSE MAB**

货号: N261032

产品全名: GTF2H1 (3E8) 小鼠单抗

基因符号 P62; BTF2; TFB1; TFIIH.

**UNIPROT ID:** P32780

背景: Component of the core-TFIIH basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II.

抗原: Purified recombinant human GTF2H1 protein fragments expressed in E.coli.

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 3E8-D3-H8

分子量: Calculated MW: 62 kDa; Observed MW: 62 kDa

亚型: IgG1


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling


储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of GTF2H1 (3E8) in HeLa using GTF2H1 antibody.



Western blot analysis of GTF2H1 in HeLa, HeLa NE and Jurkat lysates using GTF2H1 antibody.



Immunoprecipitation analysis of GTF2H1 (3E8) in HeLa lysates using GTF2H1 antibody.