

GTPASE HRAS RABBIT PAB

货号: N225374

产品全名: GTPase HRAS 兔多抗

基因符号 NS6; CMNS; NCMS; ALPS4; N-ras; NRAS1

UNIPROT ID: P01111

背景: Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Miscellaneous Mutations which change AA 12, 13 or 61 activate the potential of Ras to transform cultured cells and are implicated in a variety of human tumors.

抗原: The antiserum was produced against synthesized peptide derived from human RASH/RASK. AA range:1-50

经过测试的应用: WB, IHC-P, ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 21 kDa; Observed MW: 21 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of N/H/KRas in T47D, 3T3, COS7 and HeLa lysates using N/H/KRas antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
