

GTPASE NRAS RABBIT MAB

货号: N262816

产品全名: GTPase NRAS 兔单克隆抗体

基因符号 NS6; CMNS; NCMS; ALPS4; N-ras; NRAS1

UNIPROT ID: P01111

背景: Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Miscellaneous Mutations which change AA 12, 13 or 61 activate the potential of Ras to transform cultured cells and are implicated in a variety of human tumors.

抗原: A synthetic peptide of human Ras

经过测试的应用: WB, IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-8E1

分子量: Calculated MW: 21 kDa; Observed MW: 21 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Ras in Hela, 3T3, C6, rat Brain, MCF-7 lysates using GTPase NRAS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
