

## **GUCY2C (ICTCAR-CRC) MAB**

货号: 28581

产品全名: GUCY2C (ICTCAR-CRC) 单抗

基因符号 DIAR6, GC-C, GUC2C, MECIL, MUCIL, STAR

描述: GUCY2C (ICTCAR-CRC) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA, Flow Cyt

推荐稀释比: ELISA 1:5000-10000, Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 GUCY2C

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 $\mu$ g/mL Anti-GUCY2C (ICTCAR-CRC) mAb (28581) on Expi293 cells transfected with Human GUCY2C protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).