

HGS RABBIT PAB

货号: N225173

产品全名: HGS 兔多抗

基因符号 HGS; HRS; Hepatocyte growth factor-regulated tyrosine kinase substrate; Hrs; Protein pp110

UNIPROT ID: O14964

背景: Involved in intracellular signal transduction mediated by cytokines and growth factors.

抗原: The antiserum was produced against synthesized peptide derived from the Internal region of human HGS. AA range:511-560

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 86 kDa; Observed MW: 86 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HGS in mouse eye, mouse kidney, mouse testis lysates using HGS antibody.